

Phenol, Liquefied  
AR®  
(Contains no Preservative)



Material No.: 0025-06  
Batch No.: 0000024056  
Manufactured Date: 2011/02/24  
Retest Date: 2016/02/23

## Certificate of Analysis

Test	Specification	Result
Congealing Point	10.5 - 14.0 °C	12.7
Nonvolatile Matter	<= 0.05 %	< 0.02
Reaction	Passes Test	PT
Solubility (5 in 75)	Passes Test	PT

For Laboratory, Research or Manufacturing Use  
Appearance (clear, colorless to yellow solution)

Country of Origin: US  
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004  
Paris, KY 9001:2008  
Mexico City, Mexico 9001:2008  
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003  
Gliwice, Poland 9001:2008, 17025:2005  
Selangor, Malaysia 9001:2008  
Dehradun, India, 9001:2008, 14001:2004, 13485:2003  
Mumbai, India, 9001:2008, 17025:2005  
Panoli, India 9001:2008

Richard M Siberski  
Global Director of Quality Assurance

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600  
Avantor™ Performance Materials Inc.

3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610