

Phenol, Loose Crystals
AR® (ACS)



Material No.: 0028-10
Batch No.: 0000017492
Manufactured Date: 2012/07/27
Retest Date: 2017/07/26

Certificate of Analysis

Meets ACS Reagent Chemical Requirements,

Test	Specification	Result
ACS - Assay (C ₆ H ₅ OH)	>= 99.0 %	99.6
ACS - Appearance (white deliquescent crystals)	Passes Test	PT
ACS - Clarity of Solution	Passes Test	PT
ACS - Preservative (H ₃ PO ₂), For Information Only	%	0.05
ACS - Residue after Evaporation	<= 0.05 %	< 0.05
ACS - Water (H ₂ O)(by Karl Fischer titrn)	<= 0.5 %	0.4
ACS - Freezing Point (dry basis)	>= 40.5 °C	40.5

For Laboratory, Research, or Manufacturing Use

Country of Origin: US
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003
Gliwice, Poland 9001:2008, 17025:2005
Selangor, Malaysia 9001:2008
Dehradun, India, 9001:2008, 14001:2004, 13485:2003
Mumbai, India, 9001:2008, 17025:2005
Panoli, India 9001:2008

Richard M Siberski
Global Director of Quality Assurance

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600
Avantor™ Performance Materials Inc.

3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610