

Ethylene Dichloride
BAKER



Material No.: 9301-03
Batch No.: 0000097725
Manufactured Date: 2013/08/20
Retest Date: 2018/08/19

Certificate of Analysis

Test	Specification	Result
Specific Gravity at 25°/25°C	1.245 - 1.255	1.250
Acidity	<= 10 ppm	< 10
Free Halogens	Passes Test	PT
Lead (Pb)	<= 1 mg/kg	< 1
Nonvolatile Residue	<= 0.002 %	< 0.001
Water (H ₂ O)(by Karl Fischer titrn)	<= 0.03 %	0.02
Distilling Range(Lower)	82 - 85 °C	83
Distilling Range(Lower)	82 - 85 °C	83

For Laboratory, Research or Manufacturing Use

Country of Origin: US
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003
Gliwice, Poland 9001:2008, 17025:2005
Selangor, Malaysia 9001:2008
Dehradun, India, 9001:2008, 14001:2004, 13485:2003
Mumbai, India, 9001:2008, 17025:2005
Panoli, India 9001:2008

Richard M Siberski
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600
Avantor™ Performance Materials Inc.

3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610