

Ethylene Dichloride  
BAKER



Material No.: 9301-01  
Batch No.: 0000137228  
Manufactured Date: 2016/01/27  
Retest Date: 2021/01/25

## Certificate of Analysis

Test	Specification	Result
Specific Gravity at 25°/25°C	1.245 – 1.255	1.250
Acidity	<= 10 ppm	< 10
Free Halogens	Passes Test	PT
Lead (Pb)	<= 1 mg/kg	< 1
Nonvolatile Residue	<= 0.002 %	< 0.001
Water (H <sub>2</sub> O)(by Karl Fischer titrn)	<= 0.03 %	< 0.01
Distilling Range(Lower)	82 – 85 °C	83
Distilling Range(Lower)	82 – 85 °C	83

For Laboratory, Research or Manufacturing Use

Country of Origin: US  
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004, FSSC 22000  
Paris, KY 9001:2008  
Mexico City, Mexico 9001:2008  
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003  
Gliwice, Poland 9001:2008  
Selangor, Malaysia 9001:2008  
Dehradun, India, 9001:2008, 14001:2004, 13485:2003  
Mumbai, India, 9001:2008, 17025:2005  
Panoli, India 9001:2008

  
Jamie Ethier  
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600  
Avantor™ Performance Materials Inc.

3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610