

Hydrogen Peroxide 3% Solution, Sterile
BAKER
Made with USP Purified Water



Material No.: P006-02
Batch No.: 0000178489
Manufactured Date: 2017/08/23
Expiration Date: 2019/07/31

Certificate of Analysis

GMP Manufactured Product

Test	Specification	Result
USP – Identification	Passes Test	PT
USP – Acidity	Passes Test	PT
USP – Nonvolatile Residue (USP),mg.	<= 30	<6
USP – Barium (Ba)	Passes Test	PT
USP – Heavy Metals (as Pb)	<= 5 ppm	< 1
USP – Limit of Preservative (in mg)	<= 50	<50
USP – Assay (g/100mL)	2.5 – 3.5	3.3

CAUTION: For Manufacturing, processing or repackaging
The above referenced lot has passed the current USP Test for Sterility.

Country of Origin: US
Packaging Site: Paris Mfg Ctr & DC
Manufacturer: P0128001
Manufacturer source batch: 0000178489



Phillipsburg, NJ 9001:2008, 14001:2004, FSSC 22000
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003
Gliwice, Poland 9001:2008, 13485:2012
Selangor, Malaysia 9001:2008
Dehradun, India, 9001:2008, 14001:2004, 13485:2003
Mumbai, India, 9001:2008
Panoli, India 9001:2008

James Ethier
Jamie Ethier
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600
Avantor Performance Materials, LLC.

3477 Corporate Parkway. Center Valley, PA 18034. U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610