

Ethylene Dichloride
BAKER



Material No.: 9301-01
Batch No.: 0000218491
Manufactured Date: 2018/06/28
Retest Date: 2023/06/27
Revision No: 1

Certificate of Analysis

| Test | Specification | Result |
|---|---------------|---------|
| Specific Gravity at 25°/25°C | 1.245 – 1.255 | 1.250 |
| Acidity | <= 10 ppm | < 10 |
| Free Halogens | Passes Test | PT |
| Lead (Pb) | <= 1 mg/kg | < 1 |
| Nonvolatile Residue | <= 0.002 % | < 0.001 |
| Water (H ₂ O)(by Karl Fischer titrn) | <= 0.03 % | 0.01 |
| Distilling Range(Lower) | 82 – 85 °C | 83 |
| Distilling Range(Lower) | 82 – 85 °C | 83 |

For Laboratory, Research or Manufacturing Use

Country of Origin: US
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2015, FSC22000
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Gliwice, Poland 9001:2015, 13485:2012
Selangor, Malaysia 9001:2008
Dehradun, India, 9001:2008, 14001:2004, 13485:2003
Mumbai, India, 9001:2015, 17025:2005
Panoli, India 9001:2015


Jamie Ethier
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.386.1700
Avantor Performance Materials, LLC
100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone: 610.386.1700