

Ethylene Dichloride
BAKER



Material No.: 9301-01
Batch No.: 0000191947
Manufactured Date: 2017/03/13
Retest Date: 2022/03/12
Revision No: 1

Certificate of Analysis

Test	Specification	Result
Specific Gravity at 25°/25°C	1.245 – 1.255	1.249
Acidity	<= 10 ppm	1
Free Halogens	Passes Test	PT
Lead (Pb)	<= 1 mg/kg	< 1
Nonvolatile Residue	<= 0.002 %	< 0.001
Water (H ₂ O)(by Karl Fischer titrn)	<= 0.03 %	< 0.01
Distilling Range(Lower)	82 – 85 °C	83
Distilling Range(Lower)	82 – 85 °C	83

For Laboratory, Research or Manufacturing Use

Country of Origin: US
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004, FSSC 22000
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003
Gliwice, Poland 9001:2008, 13485:2012
Selangor, Malaysia 9001:2008
Dehradun, India, 9001:2008, 14001:2004, 13485:2003
Mumbai, India, 9001:2008
Panoli, India 9001:2008

James Ethier
Jamie Ethier
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600
Avantor Performance Materials, LLC.

3477 Corporate Parkway. Center Valley, PA 18034. U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610