

Ethylene Dichloride  
F.C.C.



Material No.: 9301-01  
Batch No.: 0000082593  
Manufactured Date: 2013/08/20  
Retest Date: 2018/08/19

## Certificate of Analysis

Meets F.C.C. Requirements,  
Food GMP Manufactured Product

Test	Specification	Result
Specific Gravity at 25°/25°C	1.245 - 1.255	1.250
Acidity	<= 10 ppm	< 10
Free Halogens	Passes Test	PT
Lead (Pb)	<= 1 mg/kg	< 1
Nonvolatile Residue	<= 0.002 %	< 0.001
Water (H <sub>2</sub> O)(by Karl Fischer titrn)	<= 0.03 %	0.02
Distilling Range(Lower)	82 - 85 °C	83
Distilling Range(Lower)	82 - 85 °C	83

Bulk Food Chemical  
CAUTION: For Manufacturing, processing or repackaging

Country of Origin: US  
Packaging Site: Paris Mfg Ctr & DC  
Manufacturer: P0005116



Phillipsburg, NJ 9001:2008, 14001:2004  
Paris, KY 9001:2008  
Mexico City, Mexico 9001:2008  
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003  
Gliwice, Poland 9001:2008, 17025:2005  
Selangor, Malaysia 9001:2008  
Dehradun, India, 9001:2008, 14001:2004, 13485:2003  
Mumbai, India, 9001:2008, 17025:2005  
Panoli, India 9001:2008

Richard M Siberski  
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600  
Avantor™ Performance Materials Inc.

3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610