

Ethylene Dichloride  
BAKER



Material No.: 9301-03  
Batch No.: 0000170150  
Manufactured Date: 2017/03/01  
Retest Date: 2022/02/28

## Certificate of Analysis

Test	Specification	Result
Specific Gravity at 25°/25°C	1.245 - 1.255	1.250
Acidity	<= 10 ppm	< 10
Free Halogens	Passes Test	PT
Lead (Pb)	<= 1 mg/kg	< 1
Nonvolatile Residue	<= 0.002 %	< 0.001
Water (H <sub>2</sub> O)(by Karl Fischer titrn)	<= 0.03 %	< 0.01
Distilling Range(Lower)	82 - 85 °C	83
Distilling Range(Lower)	82 - 85 °C	83

For Laboratory, Research or Manufacturing Use

Country of Origin: US  
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004, FSSC 22000  
Paris, KY 9001:2008  
Mexico City, Mexico 9001:2008  
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003  
Gliwice, Poland 9001:2008, 13485:2012  
Selangor, Malaysia 9001:2008  
Dehradun, India, 9001:2008, 14001:2004, 13485:2003  
Mumbai, India, 9001:2008  
Panoli, India 9001:2008

  
Jamie Ethier  
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600  
Avantor Performance Materials, LLC.

3477 Corporate Parkway. Center Valley, PA 18034. U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610