

Ethylene Dichloride  
BAKER



Material No.: 9301-03  
Batch No.: 0000226410  
Manufactured Date: 2019/02/08  
Retest Date: 2024/02/07  
Revision No: 1

## Certificate of Analysis

Test	Specification	Result
Specific Gravity at 25°/25°C	1.245 - 1.255	1.249
Acidity	<= 10 ppm	< 10
Free Halogens	Passes Test	PT
Lead (Pb)	<= 1 mg/kg	< 1
Nonvolatile Residue	<= 0.002 %	< 0.001
Water (H <sub>2</sub> O)(by Karl Fischer titrn)	<= 0.03 %	< 0.01
Distilling Range(Lower)	82 - 85 °C	83
Distilling Range(Lower)	82 - 85 °C	83

For Laboratory, Research or Manufacturing Use

Country of Origin: US  
Packaging Site: Paris Mfg Ctr & DC

  
Jamie Ethier  
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.386.1700  
Avantor Performance Materials, LLC  
100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone: 610.386.1700