

BuffAR pH 4.01 Reference Solution  
StandARd®



Material No.: 0029-04  
Batch No.: 0000009962  
Manufactured Date: 2012/07/17  
Expiration Date: 2014/01/15

## Certificate of Analysis

Test	Specification	Result
pH at 25°C	4.00 – 4.02	4.01

For Laboratory, Research, or Manufacturing Use

Traceable to NIST

Product Information (not specifications):

Appearance (clear, colorless solution)

Country of Origin: US

Packaging Site: Paris Mfg Ctr & DC

ISO

Phillipsburg, NJ 9001.2008, 14001.2004  
Paris, KY 9001.2008  
Mexico city, Mexico 9001.2008  
Deventer, The Netherlands 9001.2008, 14001.2004  
Selangor, Malaysia 9001.2008  
Panoli, India 9001.2008  
Gliwice, Poland 9001.2008, 17025.2005

Richard M Siberski  
Global Director of Quality Assurance

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600  
Avantor™ Performance Materials Inc.

3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610