

Hydrogen Peroxide, 30%  
BAKER  
(Stabilized)  
Formerly C.P.



Material No.: 2189-01  
Batch No.: 0000244028  
Manufactured Date: 2019/10/25  
Expiration Date: 2021/04/24  
Revision No: 1

## Certificate of Analysis

Test	Specification	Result
Assay (H <sub>2</sub> O <sub>2</sub> )	29.0 - 32.0 %	31.6
Residue after Evaporation	<= 0.0300 %	< 0.0300
Nitrogen Compounds (as N)	<= 0.001 %	< 0.001
Phosphate (PO <sub>4</sub> )	<= 0.02 %	< 0.02
Sulfate (SO <sub>4</sub> )	<= 0.001 %	< 0.001
Chloride (Cl)	<= 5 ppm	< 5
ACS - Heavy Metals (as Pb)	<= 1 ppm	< 1
Trace Impurities - Iron (Fe)	<= 0.5 ppm	< 0.5

For Laboratory, Research or Manufacturing Use

Country of Origin: US  
Packaging Site: Paris Mfg Ctr & DC

  
Jamie Ethier  
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.386.1700  
Avantor Performance Materials, LLC  
100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone: 610.386.1700