

Hydrogen Peroxide, 30%
BAKER
(Stabilized)
Formerly C.P.



Material No.: 2189-01
Batch No.: 0000251272
Manufactured Date: 2020/01/31
Expiration Date: 2021/07/31
Revision No: 1

Certificate of Analysis

Test	Specification	Result
Assay (H ₂ O ₂)	29.0 - 32.0 %	31.6
Residue after Evaporation	<= 0.0300 %	< 0.0300
Nitrogen Compounds (as N)	<= 0.001 %	< 0.001
Phosphate (PO ₄)	<= 0.02 %	< 0.02
Sulfate (SO ₄)	<= 0.001 %	< 0.001
Chloride (Cl)	<= 5 ppm	< 5
ACS - Heavy Metals (as Pb)	<= 1 ppm	< 1
Trace Impurities - Iron (Fe)	<= 0.5 ppm	< 0.5

For Laboratory, Research or Manufacturing Use

Country of Origin: US
Packaging Site: Paris Mfg Ctr & DC


Jamie Ethier
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.386.1700
Avantor Performance Materials, LLC
100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone: 610.386.1700